
AH Plus® & AH Plus Jet® Root Canal Sealers

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Pesquisa Medline – BC Sealer

US National Library of Medicine
National Institutes of Health

PubMed

Create RSS Create alert Advanced

Article types
Clinical Trial
Review
Customize ...

Text availability
Abstract
Free full text
Full text

Publication dates
5 years
10 years
Custom range...

Species
Humans
Other Animals

[Clear all](#)
[Show additional filters](#)

Format: Summary ▾ **Sort by:** Most Recent ▾ **Per page:** 20 ▾ **Send to** ▾

Best matches for BC Sealer AND jsubsetd:
[Push-out bond strength of root fillings made with C-Point and BC sealer versus gutta-percha and AH Plus after the instrumentation of oval canals with the Self-Adjusting File versus WaveOne.](#)
Pawar AM et al. Int Endod J. (2016)
[Antibacterial Activity of Endodontic Sealers against Planktonic Bacteria and Bacteria in Biofilms.](#)
Kapralos V et al. J Endod. (2018)
[Effect of different endodontic sealers and time of cementation on push-out bond strength of fiber posts.](#)
Vilas-Boas DA et al. Clin Oral Investig. (2018)
[Switch to our new best match sort order](#)

Search results
Items: 1 to 20 of 47
[<< First](#) [< Prev](#) Page of 3 [Next >](#) [Last >>](#)

Pesquisa Medline - Roekoseal

US National Library of Medicine
National Institutes of Health

PubMed ▼ roekoseal AND jsubsetd

Create RSS Create alert Advanced

Article types
Clinical Trial
[Review](#)
Customize ...

Text availability
Abstract
Free full text
Full text

Publication dates
5 years
10 years
Custom range...

Species
Humans
Other Animals

[Clear all](#)
[Show additional filters](#)

Format: Summary ▼ Sort by: Most Recent ▼ Per page: 20 ▼ Send to ▼

Best matches for roekoseal AND jsubsetd:
[A comparison study of periapical repair in dogs' teeth using **RoekoSeal** and AH plus root canal sealers: a histopathological evaluation.](#)
Leonardo MR et al. J Endod. (2008)
[Antibacterial Activity of Endodontic Sealers against Planktonic Bacteria and Bacteria in Biofilms.](#)
Kapralos V et al. J Endod. (2018)
[The cytotoxicity of **RoekoSeal** and AH plus compared during different setting periods.](#)
Miletić I et al. J Endod. (2005)
[Switch to our new best match sort order](#)

Search results
Items: 1 to 20 of 53

<< First < Prev Page 1 of 3 Next > Last >>

Pesquisa Medline – AH Plus (ou TopSeal ou ThermaSeal Plus)

US National Library of Medicine
National Institutes of Health

PubMed ▼ "AH Plus" OR Topseal OR Thermaseal Plus OR Ribbon Sealer AND jsubsetd

Create RSS Create alert Advanced

Article types
Clinical Trial
Review
Customize ...

Text availability
Abstract
Free full text
Full text

Publication dates
5 years
10 years
Custom range...

Species
Humans
Other Animals

[Clear all](#)
[Show additional filters](#)

Format: Summary ▼ Sort by: Most Recent ▼ Per page: 20 ▼ Send to ▼

Best matches for "AH Plus" OR Topseal OR Thermaseal Plus OR Ribbon Sealer AND jsubsetd:
[Evaluation of Physicochemical Properties of a New Root Canal Sealer.](#)
Vertuan GC et al. J Endod. (2018)
[Evaluation of the Cytotoxicity and Biocompatibility of New Resin Epoxy-based Endodontic Sealer Containing Calcium Hydroxide.](#)
Cintra LTA et al. J Endod. (2017)
[Push-out bond strength of root fillings made with C-Point and BC sealer versus gutta-percha and AH Plus after the instrumentation of oval canals with the Self-Adjusting File versus WaveOne.](#)
Pawar AM et al. Int Endod J. (2016)

Switch to our new best match sort order

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Schäfer E, Evaluation of current root canal filling materials, ZM 93:1, 24-28, article in German.

Attin, T, Rödиг, T: The root canal sealers AH 26 and AH Plus - an overview, article in German, to be published

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
Ørstavik D, Nordahl I, Tibballs JE: Dimensional change following setting of root canal sealer materials. Dent Mater 2001; 17: 512-519.

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
McMichen FR, Pearson G, Rahbaran S, Gulabivala K: A comparative study of selected physical properties of five root-canal sealers. Int Endod J 2003; 36: 629-635.

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

Força Push-out

An In Vitro Comparison of Bond Strength of Various Obturation Materials to Root Canal Dentin Using a Push-Out Test Design

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- Wennberg A: Biological evaluation of root canal sealers using in vitro and in vivo methods. J Endod 1980; 6: 784-787.
- Tagger M, Tagger E: Subcutaneous reactions to implantation of tubes with AH26 and Grossman's sealer. Oral Surg Oral Med Oral Pathol 1986; 62: 434-440.
- Figueiredo JAP, Pesce HF, Gioso MA, Figueiredo MAZ: The histological effects of four endodontic sealers implanted in the oral mucose: submucous injection versus implant in polyethylene tubes. International Endodontic Journal 2001; 34: 377-385.
- Dayan D, Buchner A, Moscona D, Sperling I: Pigmentation of the oral mucosa after root canal filling with AH-26. A light and electron microscopic study. Clin Prev Dent 1983; 5: 25-29.
- Leonardo MR, da Silva LA, Almeida WA, Utrilla LS: Tissue response to an epoxy resin-based root canal sealer. Endod Dent Traumatol 1999; 15: 28-32.
- Berbert FL, Leonardo MR, Silva LA, Tanomaru FM, Bramante CM: Influence of root canal dressings and sealers on repair of apical periodontitis after endodontic treatment. Oral Surg Oral Med Oral Pathol Oral Radiol Endod 2002; 93: 184-189.
- Pascon EA, Leonardo MR, Safavi K, Langeland K: Tissue reaction to endodontic materials: methods, criteria, assessment, and observations. Oral Surg Oral Med Oral Pathol 1991; 72: 222-237.
- Ørstavik D, Mjör IA: Usage test of four endodontic sealers in Macaca fascicularis monkeys. Oral Surg Oral Med Oral Pathol 1992; 73: 337-344.
- Bernáth M, Szabó J: Tissue reaction initiated by different sealers. Int Endod J 2003; 36: 256-261.

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Attin, T, Rödiger, T: The root canal sealers AH 26 and AH Plus - an overview, article in German, to be published

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