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Palodent Plus Clinical Case from Dr. Dao (F)

Direct composite reconstructions using the Palodent Plus[®], Sectional Matrix System

Class II direct proximal-occlusal restorations can be a challenge, as the morphology of the crown has to be reproduced to restore normal physiological function. The Palodent Plus[®] system (DENTSPLY) makes it easier to perform this type of reconstruction to preserve the dental and periodontal structures.

Introduction

Reconstruction of the crown morphology using a direct technique after caries lesions or coronal fracture is an everyday procedure in dental practices. However, class II proximal-occlusal restorations are still challenging to perform, and they must restore the crown morphology to allow normal physiological function [1].

An interproximal contact point that is not tight will not prevent food from entering the interdental space, between the interdental papilla and the interdental contact point [2], leading to accumulation of microorganisms. The presence of these microorganisms encourages caries or periodontal problems by formation of dental plaque on the tooth surface, i.e. local periodontitis [3][4]. The role of the interdental papilla is therefore to protect the dental and periodontal structures; so preserving the anatomy of the tooth, with tight contact points, is the best way of preventing caries and the best form of protection for the surrounding periodontium.

Many matrix systems are available to help dental surgeons recreate the proximal-occlusal anatomy of carious or fractured teeth, and they include the devices needed to separate the teeth [5][6]. The Palodent Plus[®] system (DENTSPLY) is one of these. This system has been redesigned and improved to allow practitioners to construct proximal contact points as simply as possible.

The Palodent Plus[®], Sectional Matrix System

The components of the Palodent Plus[®] system (DENTSPLY) are:


- Wedges, with or without protection for the adjacent tooth
- Sectional matrices with tabs
- Nickel-titanium retaining rings
- Instruments: placement forceps and pin tweezers.

A number of improvements have made the Palodent Plus[®] system very easy to use in everyday practice.

First, all the wedges and matrices have tabs or holes that make it very easy to position and place them using the pin tweezers.

Fig 1: Removing the wedge using Palodent Plus® pin tweezers.

Figs 2 and 3: Removing the matrix using Palodent Plus® pin tweezers.

There are three holes in each matrix to allow it to be removed after coronal reconstruction by sliding either occlusally or laterally. This is a very effective and practical system, as after the nickel-titanium ring has been removed, the two teeth tighten against the matrix, and it is often easier to remove the matrix by sliding it laterally.

The matrices come in different sizes, and they have also been redesigned to produce more anatomically-shaped reconstructions. They are concave, with a more pronounced edge for reconstructing the marginal ridge. The Palodent Plus® system also includes a matrix with a convex profile in its apical section, allowing more anatomical reconstruction when there is a deep defect extending down to the cementoenamel junction, which is very convenient for deep reconstructions.

The main improvement in this system is the Palodent Plus® ring. Its apical section is V-shaped, and made of glass fibre-reinforced plastic. This shape makes it much easier to place the wedges.

With older systems it was not possible to place the ring exactly where it was needed. It could be placed above the wedge, which often made it unstable (see figures below). It could also be in contact with the tooth requiring reconstruction or the adjacent tooth (the latter not allowing the best contact between matrix and tooth), but it had less effect as the forces exerted were partly absorbed by the wedge.

Fig 4: Separator ring applied over the wedge.

Fig 5: Separator ring applied distally to the wedge.

Fig 6: Separator ring applied mesially to the wedge.

So the V-shaped Palodent Plus® ring is easy to place, leaving room for the wedge, which may be inserted before or after ring placement. Both ends of the V on the ring fit precisely against the matrix on the tooth to be reconstructed. The forces they exert are thus perfectly distributed to separate the teeth to compensate for the thickness of the matrix or matrices.

Figs 7 and 8: Palodent Plus® ring in position with the Palodent Plus® wedge (vestibular and palatal view).

Clinical case studies

Clinical case study 1

The patient had a first molar with a mesial coronal fracture. The molar had an occlusal reconstruction with an infiltrated amalgam.

This was an indication for treating the mesial carious lesion and replacing the amalgam. A rubber dam was placed, the amalgam removed and the carious lesion cleaned. The Palodent Plus® matrix was then put in place and a wedge inserted in the interdental space, allowing optimum adaptation of the matrix to the tooth to be reconstructed, as shown in Fig 4. The nickel-titanium ring was then put in position, making use of the elasticity of the periodontal ligament to separate the adjacent teeth, thereby compensating for the matrix thickness. A direct reconstruction was performed using the Ceram•X® duo (DENTSPLY) composite system; the mesial wall was reconstructed to form a class I cavity. SDR® (DENTSPLY) was applied to the floor of the cavity preparation, then covered with a layer of dentine-colored composite followed by a layer of enamel-colored composite. The rubber dam was then removed, occlusion was checked and the reconstruction was polished.

Clinical case study 1

- 1 Preoperative situation.
- 2 Cleaning the cavity.
- 3 Placing the matrix.
- 4 Placing the wedge.
- 5 Placing the nickel-titanium ring.
- 6 Constructing the mesial wall (Ceram•X® duo).
- 7 Direct reconstruction with composite (Ceram•X® duo).
- 8 Checking occlusion and polishing

Clinical case study 2

This patient had two mandibular molars with proximal caries, and both teeth had infiltrated occlusal composites.

It was decided to treat both teeth at the same time. The caries cavities were cleaned thoroughly down to healthy and/or reactive dentine. Two matrices were then applied and held in place by an apical wedge. A nickel-titanium ring was then applied. Its role is even more important in this clinical situation, as the force that it will exert has to compensate for two matrices rather than one. The reconstruction was performed using Ceram•X[®] duo (DENTSPLY) and SDR[®] (DENTSPLY) flowable composite on the floor of the cavity preparation. It can be seen that the marginal ridges and proximal walls of both these cavities have been reconstructed with an anatomical shape, and that the nickel-titanium ring has exactly compensated for the thickness of both matrices, as is shown by the reconstructed contact point.

Clinical case study 2

- 1 Preoperative situation.
- 2 Cleaning the cavities.
- 3 Placing the matrices and wedge.
- 4 Placing the nickel-titanium ring.
- 5 Constructing the proximal walls (Ceram•X[®] duo).
- 6 Direct reconstruction with composite (Ceram•X[®] duo).
- 7 Checking occlusion and polishing.

Clinical case study 3


This situation was similar to the previous one, involving a patient with two carious lesions in two maxillary premolars.

It was again decided to treat both teeth at the same time. The caries lesions were cleaned and two matrices were placed, followed by a wedge and a nickel-titanium ring. The reconstruction was performed using Ceram•X[®] duo (DENTSPLY) composites. Proximal walls were reconstructed to transform both class II cavities into class I cavities. The occlusal anatomy was then reconstructed. Occlusion was checked and the reconstructions were polished.

Again, this case demonstrates preservation of the anatomy of the proximal walls and marginal ridges, with a tight contact point despite the simultaneous use of two matrices.

Clinical case study 3

- 1 Preoperative situation.
- 2 Placing the matrices.
- 3 Placing the wedge.
- 4 Placing the nickel-titanium ring.
- 5 Constructing the proximal walls (Ceram•X[®] duo).
- 6 Direct reconstruction with composite (Ceram•X[®] duo).
- 7 Direct reconstruction completed.
- 8 Checking occlusion and polishing.

Conclusion

When treating caries lesions, reconstructing the contact points between teeth is crucial for preserving dental and periodontal health. Even though these reconstructions are an everyday procedure for dental surgeons, they are still difficult to perform as many factors have to be taken into account, i.e. the anatomical morphology of the proximal walls and marginal ridges to allow dental floss to pass through, no excess composite in the reconstruction, a good matrix seal for an optimal reconstruction, and a tight contact point.

The Palodent Plus® (DENTSPLY) system makes it possible to perform class II coronal reconstructions more easily, satisfying all the requirements of these reconstructions; the anatomical matrices help to preserve tooth morphology, while the wedges make it possible to obtain a tight seal in the apical region, avoiding any excess composite. And finally, the nickel-titanium rings fulfill a number of functions, such as sealing the cavity, holding the wedge in place and making it possible to obtain a tight contact point to prevent food impaction, even when two matrices have been used simultaneously.

References

- [1] Loomans BA, Opdam NJ, Roeters FJ, Bronkhorst EM, Burgersdijk RC, Dörfer CE. A randomized clinical trial on proximal contacts of posterior composites. *Journal of Dentistry* 2006;34:292–7.
- [2] Ash MM. *Wheeler's dental anatomy physiology and occlusion*. 7th ed. Philadelphia: W.B. Saunders; 1993.
- [3] Abrahams H, Kopczyk RA. Gingival sequela from retained piece of dental floss. *Journal of the American Dental Association* 1983;106:57–8.
- [4] Hancock EB, Mayo CV, Schwab RR, Wirthlin MR. Influence of interdental contacts on periodontal status. *Journal of Periodontology* 1980;51:445–9.
- [5] Van der Vyver PJ. Posterior composite resin restorations. Part 3. Matrix systems. *Journal of the South African Dental Association* 2002;57:221–6.
- [6] ROTHE F. Coll. La reconstruction du point de contact par les matériaux plastiques. *Evaluation des mouvements, mise en œuvre*. *Odontologia*, 1985 ; 6, (2), 9-33

Original french text : published in « CLINIC » (October 2012)

Dr Dao gave authorization (email Florian Piquet, PM in France) that we can use pictures in our flyers and publish the articles internationally

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