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Posterior Class II composite restorations: The Palodent[®] Plus system

Víctor Alonso de la Peña, DMD

Associate Professor, Integrated Adult Odontology, Faculty of Medicine and Odontology, University of Santiago de Compostela, Spain.

When placing proximal restorations, anatomically shaped contact areas are important to avoid food impaction and subsequent gingival inflammation. Recent publications have confirmed that sectional matrix systems are more suitable than circumferential matrix systems for restoring class II lesions with composite in the posterior region.^{1.2}

This report describes the restoration of two upper right premolars using the sectional Palodent[®] Plus matrix system.

The preoperative radiograph showed carious lesions distally on tooth 14 and mesially and distally on tooth 15. For working comfort, if at all possible, the clamp of the rubber dam is always attached one or two teeth distally of the site we are working on so that the wings of the clamp do not interfere with the placement of the wedges and matrices (Figs. 1 and 2).

Before starting to remove the decayed tissue, we inserted protective wedges to help avoid any damage to the adjacent tooth wall. Even if the protected tooth already has a carious lesion, this is still advisable to prevent damage to sound enamel. We made sure that the teeth were spaced further apart than the thickness of the matrix (Fig. 3).

We then removed occlusal enamel with a turbine bur and any carious dentin with a round bur in an angled handpiece. The size of the bur should always be commensurate with the size of the cavity. We used a caries detector to make sure only diseased tissue is removed. After the preparation, we removed the wedge and inspected the extent of the lesions in the adjacent tooth (Figs. 4 to 7).

The matrix must be inserted without pressure, as it is very thin and can easily bend or crumple when subjected to pushing movements. We inserted the wedge where the metal shield was removed. It should always be located below the cavity floor. We then

centered the sectional matrix to be of approximately identical height buccally and lingually. We also ensured that the matrix fits to the contour of the tooth closely up to the floor of the cavity, to prevent subsequent irregularities. The ring of the matrix has a special recess for the wedges. We placed the recess such as to join the matrix at the tooth walls (Figs. 8 and 9).

Then we applied X-Flow[®] at the bottom of the cavity, covering the entire cavity floor

(Fig. 10). This facilitated the subsequent application of increments of Ceram.X[®] mono M1 posterior composite. The composite inserted into the cavity does not require condensing, as the matrix serves as the proximal wall. The formation of bubbles and puts must be avoided. In the posterior region, dentin shades are never used even in deep cavities (Fig. 11).

The proximal surfaces were polished with disks whose working sides face the head of the angled handpiece. Dental floss guided along the walls, if it frays, it indicates irregularities that must be eliminated (Fig. 12).

Tooth 15 was restored in a subsequent session.

Wedging is imperative before accessing the cavity. The rings are placed with a mesial orientation, as the arch of the clamp can interfere with their placement. The recess for the wedges provide stability. Care must be taken to ensure that the matrix adheres to the gingival edge of the wall. Any debris that may have been entrapped in the sulcus must be removed (Figs. 13 to 16).

Considerations

- The spacing between the teeth must be greater than the thickness of the matrix. A wedge is efficient when pressure is needed. It is therefore preferable to use the widest possible wedges and to insert them with a mosquito clamp.
- The most critical part of the clinical sequence is to introduce the matrix without deforming it; it is very thin and may be distort by pressure. The teeth must always be sufficiently separated. If the matrix cannot be inserted properly, we have the following options:
 - Inserting a wide wooden wedge and waiting 90 seconds (Figs. 17 to 20).
 - If interference is located close to the occlusal surface, it might be pried loose with a metal instrument.

- If the wing of the clamp interferes with the positioning of the wedge or matrix, use a smaller clamp without wings or place it on a tooth further back, if possible (Fig. 21).
- If the ring encroaches on the cavity, to achieve adequate mesiodistal width, a circumferential matrix system has to be used.

References

Saber MH, El-Badrawy W, Loomans BA, Ahmed DR, Dörfer CE, El Zohairy A. Creating tight proximal contacts for MOD resin composite restorations. Oper Dent. 2011; 36:304-310.
Saber MH, Loomans BA, El Zohairy A, Dörfer CE, El-Badrawy W. Evaluation of proximal contact tightness of Class II resin composite restorations. Oper Dent. 2010; 35:37-43.

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