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Class II posterior preparations: New tools

Dr. Juan José Rodríguez Ortega

Professor of Biomaterials and Dental Equipment, Faculty of Dentistry, University of Valladolid, Spain

A 55-year-old patient with no relevant medical history presented with caries (MOD) in tooth 26. We decided to prepare the tooth using the new Palodent® Plus matrices and to place a restoration with SDR® – Smart Dentin Replacement – and occlusally CeramX® duo composite.

The tooth was prepared minimally invasive, removing the demineralised enamel and infected dentin. Care was taken not to remove affected dentin, i.e. dentin that was still firm but easily accessed with a probe. Standard cutters were used to prepare the enamel, and a round bur in an angled handpiece was used to prepare the dentin (Figure 1).

The enamel was etched with phosphoric acid 37% for 5 seconds (Figure 2). The site was rinsed and dried. XP Bond® adhesive was applied to the enamel and dentine, allowed to stay for 30 seconds and light-cured (Figure 3).

Subsequently SDR® was placed (Figure 4). Application is simple; the material is provided in convenient compulas and can be bulk-applied in thick increments (4 mm).

The new Palodent® Plus system with its wedges and rings achieved a good fit. No air bubbles or voids occurred and no imperfections were created. The time savings that SDR® provides allow the dentist to concentrate on other critical aspects of the restoration, such as the formation of labial, lingual and occlusal apertures or the design of the marginal ridge.

No excess material or voids formed between the restoration and the tooth hard tissue. The application of the SDR® base material and simultaneous restoration of the marginal ridges were carried out in just 15 minutes.

Finally, the occlusal layer was built up (Figures 5 through 7). Cusps were formed using a build-up technique, with each layer being light-cured for 3 seconds. This technique not only drastically reduces the polymerisation stress and the formation of »white lines« but also considerably shortens the working time. After the final layer, the entire restoration was again light-cured for 20 seconds. Staining material was applied to

highlight of the central fossa and primary sulci and once again light-cured for 20 seconds.

Finishing and polishing were performed using the Enhance® tools and pastes (Figures 8 and 9). The final result after occlusal adjustment and final polishing is shown in Figure 10.

Materials and methods

We used Palodent® Plus matrices and SDR® base composite (both DENTSPLY), Ceram®X duo composite (DENTSPLY), and Tetric Color dye (Ivoclar).

Picture 1

Picture 2

Picture 3

Picture 4

Picture 5

Picture 6

Picture 7

Picture 8

Picture 9

Picture 10

Dr. Juan José Rodríguez Ortega

Original published Spanish case : See below (published in Dentpsly Noticias, September2012)

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