

DENTSPLY

FIRST IN DENTISTRY



TRIED&TRUE MEETS
NEW&IMPROVED.

Palodent[®] **Plus**
Sectional Matrix System


Philippe ARDILLIER - Dr Jérôme DAO

October 2011

Introkit



1st clinical case



Preparation of cavities



inlay



Inlay put in place
with Calibra




2 matrix



matrix are maintained with plastic corners



you seal up each cavity one by one by blocking the matrix with the ring





We used SDR



We finished with Ceram mono



Palodent 'Plus

Philippe ARDILLIER - Jérôme DAO

Final result




2nd clinical case



TRIED&TRUE MEETS
NEW&IMPROVED.

Palodent[®] **Plus**⁺⁺
Sectional Matrix System

Old composite



we maintain the matrix with plastic corners and we block matrix with the ring





Palodent 'Plus




Palodent 'Plus

Philippe ARDILLIER - Jérôme DAO



We used Ceram X mono

the pictures show the ease of use forceps with a pin



Final result



Palodent 'Plus

Philippe ARDILLIER - Jérôme DAO